



Cylinder, Graduated, 1000 mL, Certified

Številka izdelka: 2636253

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Specifikacije

Certificiranje:	Class A
Kapaciteta:	1000 mL
Material ohišja:	Glass
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	465 mm
Vrsta:	eno metričen