



Cylinder, graduated, mixing, certified, 50 mL

Številka izdelka:	2636341
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Specifikacije

Certificiranje:	Class A
Kapaciteta:	50 mL
Material ohišja:	Glass
Merilna metoda:	To Contain
Prostornina / velikost paketa:	1 each
Vrsta:	eno metričen, mešan