



Flask, volumetric, class A, certified, 1000 mL

Številka izdelka: 2636653

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (Stopper size 22).

Specifikacije

Certificiranje:	Yes
Enota:	1 each
Kapaciteta:	1000 mL
Material ohišja:	Class A Glass