



Be Right™



Gas washing bottle, 1200 mL

Številka izdelka: 2662200

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis—Hach Method 10065.

Specifikacije

Body Material:	Glass
Kapaciteta:	1200 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Prostornina / velikost paketa:	1 each