



Iron (Ferrous) Color Disc Test Kit, Model IR-18C, 0-7 mg/L

Številka izdelka:

2667200

EUR Cena:

Kontakt

Odprema v 2 tednih

Nearno



The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons. Range: 0-7 mg/L. Contains 100 tests.

Simple, fast, color disc method

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Specifikacije

Ime metode:	Color disc/ 1,10 Phenanthroline
Merilno območje:	High Single Range
Model:	IR-18C
Parameter:	Iron, medium range - As Fe ²⁺ (ferrous iron)
Platforma:	Color Disc
Smallest increment:	0.2
Število testov:	100
Test Kit model:	IR-18C
Vrsta:	eno metričen
Vrsta primera:	B