



## Test kit, hypochlorite, model CN-HRDT

Številka izdelka: 2687100

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

**Nearno**



Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl<sub>2</sub>: 50 - 150 g/L as Cl<sub>2</sub>. Contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities**

**Verify concentration of hypochlorite solutions at time of delivery**

**Kit contains all required reagents, solutions, and glassware**

### Specifikacije

Ime metode:	Digital Titrator/ Thiosulfate
Merilno območje:	5 - 15 % = 50-150 g/L Cl <sub>2</sub>
Parameter:	Hypochlorite, high range (Bleach) - As Cl <sub>2</sub>
Platforma:	Digitalni Titrator
Smallest increment:	0.05%
Številka metode:	10100
Število testov:	100
Test Kit model:	CN-HRDT
Vrsta primera:	G