



SpecCheck Secondary Gel Standard Kit - Fluoride, 0-2.0 mg/L F

Številka izdelka: 2712500

EUR Cena: Kontakt
Odprema v 2 tednih

Solid, ready for use Gel-Standards for instrument verification in the lab and in the field.

Set contains four vials filled with gels that simulate the test colour at various concentrations.

Specifikacije

Ime metode:	SPADNS
Inštrument:	All except DR/3000, DR 900 and DR/800 Series
Koncentracija:	0 - 2,00 mg/L F