



## Test strips free & total chlorine, 0-10 mg/L, individually wrapped

Številka izdelka: 2793944  
EUR Cena: Kontakt  
Odprema v 2 tednih

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Chlorine exists in water as "free" chlorine (hypochlorite) or "combined" chlorine (chloramines). The total chlorine test detects both forms. With Hach Water Quality Test Strips, technicians in the field can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Hach test strips are also used in laboratories all over the world for pre-test screening-to detect the presence of materials that might interfere with lab testing. Strips are individually wrapped for maximum portability and stability. Contains 250 tests.

**Easy to use, disposable, and inexpensive**

**A great way to obtain quick, quantitative answers in the field or in the lab**

---

### Specifikacije

Enota:	250 tests (individually wrapped)
Ime metode:	Test Strips
Komentar:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	0 - 10 mg/L Cl <sub>2</sub>
Najmanjši korak rasti:	Steps: 0, 0.5, 1, 2, 4, 10
Parameter:	Chlorine, low range - As free & total Cl <sub>2</sub>
Platforma:	Test Strip
Test Kit model:	TS