



## Plastic ampules R2A, pk/50

Številka izdelka:	2812350
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

R2A Broth Ampules; Plastic

R2A Plastic Ampules for detection of stressed heterotrophic bacteria; pk/50.

### Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

---

## Specifikacije

Material:	Plastic
Media Type:	R2A
Prostornina / velikost paketa:	50 /pk