



Quality control standard hardness, high range, 500 mL

Številka izdelka:

2833349

Nearno

EUR Cena:

Kontakt

Na voljo



Ready to use - no dilutions required. Mixed standard includes: Total Hardness: 1000 mg/L as CaCO₃, Calcium Hardness: 500 mg/L as CaCO₃. 500 mL Bottle

Specifikacije

Koncentracija:

Calcium Hardness 500 mg/L CaCO₃

Opis:

High Range Hardness, Mixed Parameter QC Standard

Platforma:

Raztopina