



Quality control standard oxygen demand , pk/16 - 10 mL Voluette® ampule

Številka izdelka: 2833510

EUR Cena: Kontakt

Odprema v 2 tednih

Ready to use - no dilutions required. Oxygen Demand Quality Control Standard. BOD - 396 mg/L. COD - 613 mg/L. TOC - 242 mg/L.

Specifikacije

Koncentracija:

COD - 617 mg/L

Opis:

Oxygen Demand, Mixed Parameter QC Standard

Platforma:

Raztopina