



Quality control standard metals for drinking water, low range, 500 mL

Številka izdelka: 2833749

EUR Cena: Kontakt

Na voljo

Ready to use - no dilutions required. Mixed standard includes: Copper: 1 mg/l Cu, Iron: 0.3 mg/L Fe, Manganese: 0.1 mg/L Mn. 500 mL Bottle.

Specifikacije

Koncentracija:

Mangan: 0,1 mg/L Mn

Opis:

Metals Drinking Water (Low Range) Mixed Parameter QC Standard

Platforma:

Raztopina