



## Fluoride Standard Solution, 1 mg/L, 1 L

Številka izdelka:	29153
EUR Cena:	Kontakt
Odprema v 2 tednih	

### Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Fluoride, 1.00 mg/L  $\pm 0.02$  mg/L F<sup>-</sup> (NIST). 1 L bottle

#### Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifikacije

Koncentracija:	1,0 mg/L F <sup>-</sup>
Opis:	Fluoride Standard Solution
Package settings:	Bottle
Parameter:	Fluoride
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	1 L