



Cylinder, graduated, double metric, 100 mL

Številka izdelka: 2928302

EUR Cena: Kontakt

Odprema v 2 tednih

These glass pyrex@ graduated cylinders feature two sets of graduations for reading either up or down Calibrated "to deliver"

Specifikacije

Certificiranje:	Class A
Kapaciteta:	100 mL
Material ohišja:	Glass
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	254 mm
Vrsta:	Double Metric