



Cylinder, graduated, double metric, 250 mL

Številka izdelka: 2928303
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

These glass pyrex graduated cylinders feature two sets of graduations for reading either up or down Calibrated "to deliver"

Specifikacije

Certificiranje:	Class A
Kapaciteta:	250 mL
Material ohišja:	Glass
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	340 mm
Vrsta:	Double Metric