



Xenosep consumables kit for EPA method 1664A testing

Številka izdelka: 2943800
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Order the Xenosep Consumables Kit for a quick and easy way to replace your most used components.

Consumables kit includes: 24 Xenosep SPE Filters, 47mm diameter, 24 Transfer Pipets, 15 mL capacity per pipet, 24 Inline Columns with Sodium Sulfate, 10g per column.

Specifikacije

Prostornina / velikost paketa: 1 each