



Fixed Volume, Micropipette, 10 µL

Številka izdelka: 2946002

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and colour-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikacije

Kapaciteta:	10 µL
Natančnost:	1,0 %
Prostornina / velikost paketa:	1 each