



SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 1 L

Številka izdelka:

2947553

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F⁻. EPA compliant. 1 Liter bottle.

Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

Specifikacije

EPA compliant:	Da
Ime metode:	SPADNS 2
Merilno območje:	0.02 - 2.00 mg/L F ⁻
Parameter:	Fluoride
Platforma:	Raztopina
Številka metode:	10225
Število testov:	250