



## Conductivity Standard Solution, 10000 µS/cm, NaCl, 50 mL

Številka izdelka:	2972226
EUR Cena:	Kontakt
Odprema v 2 tednih	

### Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 10000 µS/cm  $\pm$ 100 µS/cm at 25°C, 5.556 g/L  $\pm$ 0.03 g/L as Sodium Chloride (NaCl), (1 µS/cm  $\sim$  0.5 mg/L TDS), 50 mL pour-top bottle.

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifikacije

Formula:	NaCl
Koncentracija:	10000 µS/cm at 25°C $\pm$ 100 µS/cm at 25°C
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parameter:	Conductivity
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	50 mL
Sledljivo po NIST-u:	Da