



Conductivity Standard Solution, 12890 µS/cm, KCl, 50 mL

Številka izdelka: 2974426

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Conductivity standard solutions for laboratory and field analyses

Conductivity Standard Solution, 12890 µS/cm ±10% at 25°C, 7455 mg/L ±10% as Potassium Chloride (KCl) (0.1M). 50 mL pour-top bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Formula:	KCl
Koncentracija:	0.1M KCl; 7455 mg/L (KCl)
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parameter:	Conductivity
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	50 mL
Sledljivo po NIST-u:	Da