



Natural Water Standard Solution, 3000 ppm Total Dissolved Solids (TDS), 50 mL

Številka izdelka: 2974826

EUR Cena: Kontakt

Odprema v 2 tednih

Conductivity standard solutions for laboratory and field analyses

Natural Water Standard Solution. 40% sodium sulfate/40% sodium bicarbonate/20% sodium chloride. TDS: 3000 ppm \pm 30 ppm, Conductivity: 3960 μ S/cm \pm 39.6 μ S/cm at 25°C. 50 mL pour-top bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Formula:	Na ₂ SO ₄ , Na ₂ CO ₃ , NaCl
Koncentracija:	3960 μ S/cm at 25°C \pm 39.6 μ S/cm at 25°C
Opis:	TDS/Conductivity calibration solution
Package settings:	Bottle
Parameter:	TDS
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	50 mL
Sledljivo po NIST-u:	Da