



pH paper, 1.0 - 11.0 pH range

Številka izdelka:	39133
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

Specifikacije

Merilno območje:	1.0 - 11.0 pH units
Prostornina / velikost paketa:	5 4.6-meter (15') rolls per package