



Sodium/Potassium Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Številka izdelka:

4451569

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Ionic Strength Adjustor (ISA) for Sodium/Potassium determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Sodium/Potassium level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specifikacije

Koncentracija:	ISA
Package settings:	Powder Pillows
Parameter:	Sodium
Prostornina / velikost paketa:	100 /pk