



MultiTest kit case (218 x 246 x 437 mm), blue polypropylene

Številka izdelka:	4660700
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Specifikacije

Dimenzija (v x š x g):	218 mm x 246 mm x 437 mm
Višina:	8.5 in
Vrsta primera:	Q