



MultiTest kit case (295 x 330 x 576 mm), blue polypropylene

Številka izdelka: 4660800
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Specifikacije

Dimenzija (v x š x g):	295 mm x 330 mm x 576 mm
Višina:	11 in
Vrsta primera:	J