



Beaker, griffin, low form, glass, 400 mL

Številka izdelka: 50048
EUR Cena: Kontakt
Odprema v 2 tednih

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specifikacije

Enota:	1 each
Kapaciteta:	400 mL
Material:	Glass, Griffin, Low Form
Material ohišja:	Glass