



Beaker, Griffin, Low Form, Glass, 600 mL

Številka izdelka:	50052
EUR Cena:	Kontakt
Odprema v 2 tednih	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Enota:	1 each
Kapaciteta:	600 mL
Material:	Glass, Griffin, Low Form
Material ohišja:	Glass