



Beaker, griffin, low form, glass, 50 mL, 12/pk

Številka izdelka: 50071H

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specifikacije

Enota:	12 /pk
Kapaciteta:	50 mL
Material:	Glass, Griffin, Low Form
Material ohišja:	Glass