



QL filters, pore size 8 - 12 μm , diameter 70 mm, circular, 100/ pk, qualitative

Številka izdelka: 50653

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specifikacije

Diameter:	70 mm
Material:	Cellulose
Prostornina / velikost paketa:	100 /pk