



Cylinder, graduated, 5 mL $\pm$ 0.1 mL, 0.1 mL divisions (white markings)

Številka izdelka: 50837
EUR Cena: Kontakt
Na voljo

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes $\geq$ 25 mL.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	5 mL
Material ohišja:	Glass
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	125 mm
Vrsta:	eno metričen