



Cylinder, graduated, 10 mL +- 0.1 mL, 0.2 mL divisions (white markings)

Številka izdelka: 50838
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	10 mL
Material ohišja:	Glass
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	135 mm
Vrsta:	eno metričen