



Cylinder, graduated, 50 mL +/- 0.4 mL, 1.0 mL divisions (white markings)

Številka izdelka: 50841
EUR Cena: Kontakt
Odprema v 2 tednih

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	50 mL
Material ohišja:	Glass
Merilna metoda:	To Contain
Prostornina / velikost paketa:	1 each
Višina:	190 mm
Vrsta:	eno metričen