



Cylinder, graduated, 100 mL +- 0.6 mL, 1.0 mL divisions (white markings)

Številka izdelka: 50842
EUR Cena: Kontakt
Odprema v 2 tednih

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specifikacije

Certificiranje:	Class B
Kapaciteta:	100 mL
Material ohišja:	Glass
Merilna metoda:	To Deliver
Prostornina / velikost paketa:	1 each
Višina:	255 mm
Vrsta:	eno metričen