



Container, glass, 10 L (2.5 gal), with PTFE-lined cap

Številka izdelka: 6559
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

2.5 gal glass container with cap

Specifikacije

Vsebina paketa:

Bottle(6358)

Cap (6355)

PTFE Liner (6378)

Add Required parts:

For Compact or Standard 12 Bottle Portable Base

Container Support (1502)

For Refrigerated and All Weather Samplers

Container Support (1502)

Extension Tube (3527)

Vsebina paketa

Bottle(6358) Cap (6355) PTFE Liner (6378)Add Required parts:For Compact or Standard 12 Bottle Portable BaseContainer Support (1502)For Refrigerated and All Weather SamplersContainer Support (1502)Extension Tube (3527)