



Brush, buret, 50 mL

Številka izdelka:	68600
EUR Cena:	Kontakt
Odprema v 2 tednih	

Ideally shaped for quick and thorough cleaning of burets and other similarly shaped vessels. Unbleached white nylon bristles are supported by twisted galvanized steel wire. Brush portion measures 8 cm (3.1") in length. Fan style tip. 2.0cm (0.79") Dia. 91cm (35.8") L

Specifikacije

Prostornina / velikost paketa:	1 each
--------------------------------	--------