



Buffer Powder Pillows, pH 8.30, 50 mL, pk/25

Številka izdelka:

89868

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Buffer for preparing monochloramine standards and detecting visual alkalinity endpoints. pH 8.30, Sample Size: 50 mL sample. Pack of 25 Powder Pillows.

Fresh and reliable

Unit-dose powder pillows assure buffer freshness and eliminate the possibility of contamination.

Specifikacije

Koncentracija:

pH 8,3

Opis:

Buffer Powder Pillows, 25/pk

Platforma:

Praškasti reagenti