



## Recharge kit of 4 pre-filled cartridges with premounted 29552A membranes for C1100 ozone sensors with 33051-xT cap only

Številka izdelka:

29552A-CT

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Includes o-rings and cleaning tools.

---

### Specifikacije

|                      |                                          |
|----------------------|------------------------------------------|
| Inštrument:          | Orbisphere C1100 O3 sensor               |
| Parameter:           | O3                                       |
| Pogoji skladiščenja: | Store in a dry place at room temperature |