



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Številka izdelka:

2956A-AT05

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Includes o-rings and cleaning tools. Special low concentration electrolyte.

Specifikacije

Inštrument:

Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor

Parameter:

O2

Pogoji skladiščenja:

Store in a dry place at room temperature