



## **Xenomax TSS-FVS Filter 47 mm, 100/pk**

**Številka izdelka:** 2511307

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Xenomax TSS-FVS Filters are specifically designed to ensure compliance with the filter preparation requirements of Standard Methods 2540 D and Standard Methods 2540 E (Fixed and Volatile Solids Ignited at 550 °C). Independent 3rd party evaluation of Xenomax and 934-AH filters in both methods give demonstrably equivalent results.

In addition, these 1.5 µm binder-free, fast flowing, fine porosity borosilicate glass microfiber filters have been optimized for better wet strength/handleability and high temperature (550 °C) blank stability. FVS Filters are available in popular diameter sizes and ready for use in a wide variety of labware including Buchner funnels and Gooch Crucibles.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare.

---

### **Specifikacije**

|                                |                  |
|--------------------------------|------------------|
| Diameter:                      | 47 mm            |
| Material:                      | Glass Microfiber |
| Prostornina / velikost paketa: | 100 /pk          |
| Velikost por:                  | 1.5 µm           |