



CL17sc Cell Cleaning Kit

Številka izdelka: 8573100

EUR Cena: Kontakt

Na voljo

Nearno



Hach recommends cleaning the CL17sc cell on a monthly (or as needed) basis. Frequency of cleaning may depend on the application and sample qualities.

The CL17sc Cell Cleaning Kit contains 5.25N sulfuric acid (100 mL dropper bottle) and cotton swabs (10).

Specifikacije

Inštrument:	CL17sc
Merilno območje:	N/A
Parameter:	N/A
Pogoji skladiščenja:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Vsebina paketa:	Sulfuric acid (100 mL dropper bottle) and cotton swabs (10).

Vsebina paketa

Sulfuric acid (100 mL dropper bottle) and cotton swabs (10).