



CL17sc Calibration Verification Kit

Številka izdelka:	8568200
EUR Cena:	Kontakt
Na voljo	

The CL17sc Calibration Verification Kit includes two bottles of deionized water (500 mL), a chlorine standard (20 mL ampule), ampule breaker, syringe (50 mL), tubing, and fittings.

The CL17sc calibration curve is set at the factory for performance to specifications. User adjustment of the factory calibration curve is only recommended in two instances: if required by a regulatory agency for compliance reporting purposes or if a large repair of the analyzer is done.

Specifikacije

Inštrument:	CL17sc
Merilno območje:	N/A
Parameter:	N/A
Vsebina paketa:	2x 500 mL bottle of deionized water, 20 mL ampule of chlorine standard, ampule breaker, 50 mL syringe, tubing, and fittings.

Vsebina paketa

2x 500 mL bottle of deionized water, 20 mL ampule of chlorine standard, ampule breaker, 50 mL syringe, tubing, and fittings.