



CL17sc Calibration Verification Refill Kit

Številka izdelka: 8573200

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

For users that already have the syringe, tubing, fittings, and ampule breaker from the CL17sc Calibration Verification Kit, this refill kit provides the deionized water and chlorine standard needed to verify calibration. The CL17sc Calibration Verification Refill Kit contains two bottles of deionized water (500 mL) and a chlorine standard (20 mL ampule).

The CL17sc calibration curve is set at the factory for performance to specifications. User adjustment of the factory calibration curve is only recommended in two instances: if required by a regulatory agency for compliance reporting purposes or if a large repair of the analyzer is done.

Specifikacije

Inštrument:	CL17sc
Merilno območje:	N/A
Parameter:	N/A
Vsebina paketa:	2x 500 mL bottle of deionized water and 20 mL ampule of chlorine standard.

Vsebina paketa

2x 500 mL bottle of deionized water and 20 mL ampule of chlorine standard.