



## Acid solution for EZ4004 Total Alkalinity Analyser, 0.05 N HCl (25-1250 mg/L Range), 2 L

Številka izdelka:

APPC4004-02-EU **Nearno**

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan



### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifikacije

Koncentracija:	0.05 N HCl
Parameter:	Total Alkalinity
Prostornina / velikost paketa:	2 L