




## Ammonia Chlorine Solution for EZ1102 Ammonia Analyser, 500 mL

Številka izdelka: **APPC1102-01**

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifikacije

|                                |         |
|--------------------------------|---------|
| Parameter:                     | Ammonia |
| Prostornina / velikost paketa: | 500 mL  |