



## Ascorbic Acid for EZ1001 Aluminium Analyser, 100 g

Številka izdelka:	APPC1001-03B
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specifikacije

Parameter:	Aluminium
Prostornina / velikost paketa:	100 g