



## Cylinder, Graduated, Mixing, Glass, 50 mL $\pm$ 0.4 mL, 1.0 mL divisions, glass stopper #16

Številka izdelka: 104592

EUR Cena: Kontakt

Tall form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

---

### Specifikacije

|                                |                     |
|--------------------------------|---------------------|
| Certificiranje:                | Class B             |
| Kapaciteta:                    | 50 mL               |
| Material ohišja:               | Glass               |
| Merilna metoda:                | To Contain          |
| Prostornina / velikost paketa: | 1 each              |
| Vrsta:                         | eno metričen, mešan |