




**Be Right™**



## **Orthophosphate Standard Solution, 3 mg/L PO<sub>4</sub>-P, 500 mL**

**Številka izdelka:** LCW1053

EUR Cena: Kontakt

Na voljo

Standard Solution for use with the NP6000sc Phosphate Analyser.

---

### **Specifikacije**

|                                |                           |
|--------------------------------|---------------------------|
| Inštrument:                    | NP6000sc                  |
| Koncentracija:                 | 3 mg/L PO <sub>4</sub> -P |
| Parameter:                     | Orthophosphate            |
| Prostornina / velikost paketa: | 500 mL                    |