



**Be Right™**



## **Orthophosphate Standard Solution, 10 mg/L PO<sub>4</sub>-P, 500 mL**

**Številka izdelka:** LCW1055

EUR Cena: Kontakt

Na voljo

Standard Solution for use with the NP6000sc Phosphate Analyser.

---

### **Specifikacije**

Inštrument:	NP6000sc
Koncentracija:	10 mg/L PO <sub>4</sub> -P
Parameter:	Orthophosphate
Prostornina / velikost paketa:	500 mL