




Be Right™



## NitraVer 5 Nitrate Reagent Powder Pillows, 5 mL, pk/100

Številka izdelka:

1403599-EU

Nearno

EUR Cena:

Kontakt

Na voljo



Powder reagent in sealed packages for the determination of Nitrogen-Nitrate up to 50 mg/l NO<sub>3</sub>-N.

For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 100 powder pillows.

---

### Specifikacije

Ime metode:	Cadmium Reduction
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	100 /pk

---

### Obvezni pripomočki

- Zaščitne rokavice velikosti L, modre, nitril, brez pudra, 50 kosov (Item SM743L)