



Be Right™



Chlorine, free; colour disc test kit, model CN-66F (pk/100)

Številka izdelka: 223102-EU

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Measuring range 0.1 - 3.4 mg/L Cl₂

Specifikacije

Ime metode:	Color disc/DPD
Komentar:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Merilno območje:	0 - 3.4 mg/L Cl ₂
Model:	CN-66F
Parameter:	Chlorine, low range - As free Cl ₂
Platforma:	Color Disc
Smallest increment:	0.2
Standardno merilno območje:	0 - 3.4 mg/L Cl ₂
Število testov:	100
Test Kit model:	CN-66F ²
Vrsta primera:	B